

CHAMOT

重组人成纤维细胞生长因子-11亚型2 /Recombinant FGF-11 isoform 2, Human, Animal-Free

CM101-5HP

CM101-20HP

CM101-100HP

CM101-500HP

CM101-1000HP



CHAMOT

乔默®生物

Specialize In Cytokines

CONTENT

1 产品简介

2 产品使用

3 产品储存/运输

4 实验数据展示

重组人成纤维细胞生长因子-11亚型2 /Recombinant FGF-11 isoform 2, Human, Animal-Free

产品编号	CM101-5HP	CM101-20HP	CM101-100HP	CM101-500HP	CM101-1000HP
规格	5 µg	20 µg	100 µg	500 µg	1 mg

产品简介

背景描述

The Fibroblast Growth Factors (FGF) family has two group, secreted FGFs and intracellular FGFs (iFGFs). iFGFs are related to voltage-gated sodium (Nav) channels. Human FGF-11 isoform 2 is a 18 kDa cytokine with 166 amino acid residues.

别称

fibroblast growth factor 11 isoform 2, FHF-3, FHF3, FGF11

蛋白编码

Q92914

分子量

The protein has a calculated MW of 19.34 kDa. The protein migrates as 23 kDa under reducing condition (SDS-PAGE analysis).

表达系统

Escherichia coli

纯度

>98% as determined by SDS-PAGE.

生物活性

Measure by its ability to induce 3T3 cells proliferation. The ED₅₀ for this effect is <1 ng/mL.

内毒素检测

<0.1 EU per 1 µg of the protein by the LAL method.

蛋白序列

MSLSPEPQLKGIVTKLFCRQGFFYLQANPDGSIQGTPEDTSSFTHFNLIPVGLRVVTIQSAKLGHY
MAMNAEGLLYSSPHFTAECRFKECVFENYYVLYASALYRQRRSGRAWYLGDKQVMKGNR
VKKTKAAAHFLPKLLEVAMYQEPSLHSVPEASPSSPPAP with polyhistidine tag at the C-terminus.

产品形式

The protein was lyophilized from a 0.2 µm filtered solution containing 1X PBS, pH 7.4. If you have any concerns or special requirements, please confirm with us.

产品应用

Cell Culture

产品使用

1. Before opening, centrifuge at 3000 rpm for 5 mins.

2. Initial Reconstitution

- Reconstitute the Lyophilized Protein in sterile H₂O to a concentration of 100-200 µg/mL.
- Then, incubate it at room temperature for at least 20 mins to ensure sufficient dissolution.

Do Not Vortex! Vigorous shaking may impair the biological activity of the protein.

- Store at 2°C to 8°C for up to 1 week.

3. Extended Storage

- After the initial reconstitution, further dilute the reconstituted protein in a buffer containing a carrier protein or stabilizer (e.g., 0.1% BSA, 10% FBS, 5% HSA, or 5% trehalose solution). The final concentration is not less than 10 µg/ml.
- Prepare aliquots (≥20 µl).
- Store at -20°C or -80°C for 3 to 6 months.

Avoid repeated freeze-thaw cycles.

产品储存/运输

Lyophilized Protein: Store at -20°C for 1 year.

储存



Initial Reconstitution: Store at 2°C to 8°C for up to 1 week.

Extended Storage: Store at -20°C or -80°C for 3 to 6 months with a carrier protein or stabilizer.

运输

Blue Ice

实验数据展示



SDS-PAGE analysis of recombinant human FGF-

11 isoform 2

Human FGF-11 isoform 2, His Tag, E. coli induced 3T3 cell

proliferation, with the ED₅₀ at 0.84 ng/mL.

For Research Use or Further Manufacturing Only

Chamot Biotechnology(Shanghai) Co., Ltd. www.chamot-bio.com

Tel: 021-51880030 Mail: info@chamot-bio.com QQ: 864920491