

CHAMOT

**重组人白细胞介素28A /Recombinant IL-28A,
Human, Animal-Free**

CM032-5HP

CM032-20HP

CM032-100HP

CM032-500HP

CM032-1000HP



CHAMOT

乔默®生物

Specialize In Cytokines

CONTENT

1 产品简介

2 产品使用

3 产品储存/运输

4 实验数据展示

重组人白细胞介素28A /Recombinant IL-28A, Human, Animal-Free

产品编号	CM032-5HP	CM032-20HP	CM032-100HP	CM032-500HP	CM032-1000HP
规格	5 µg	20 µg	100 µg	500 µg	1 mg

产品简介

背景描述	Interleukin 28A (IL-28A) predicts a molecular mass of 22.3 kDa. IL-28 comes in two isoforms, IL-28A and IL-28B, and plays a role in immune defense against viruses, including the induction of an antiviral state by turning on Mx proteins, 2',5'-oligoadenylate synthetase as well as ISGF3G.
别称	interleukin 28A, IFN-λ2, IL28, IL28A
蛋白编码	Q8IZJ0
分子量	The protein has a calculated MW of 20.60 kDa. The protein migrates as 20 kDa under reducing condition (SDS-PAGE analysis).
表达系统	Escherichia coli
纯度	>98% as determined by SDS-PAGE.
生物活性	Measure by its ability to protect HepG2 cells infected with encephalomyocarditis (EMC) virus. The ED ₅₀ for this effect is <20 ng/mL.
内毒素检测	<0.1 EU per 1 µg of the protein by the LAL method.
蛋白序列	MPVARLHGALPDARGCHIAQFKSLSPQELQAFKRAKDALEESLLLKDCRCHSRLFPRTWDLRQLQVRERPMALAEALATLKVLEATADTDPALVDVLDQPLHTLHHILSQFRACIQPQPTAGPRTRGRLHHWLYRLQEAPKKESPGCLEASVTFNLFRLLTRDLNLCVASGDLCV with polyhistidine tag at the C-terminus.
产品形式	The protein was lyophilized from a 0.2 µm filtered solution containing 1X PBS, pH 8.0. If you have any concerns or special requirements, please confirm with us.
产品应用	Cell Culture

产品使用

1. Before opening, centrifuge at 3000 rpm for 5 mins.

2. Initial Reconstitution

- Reconstitute the Lyophilized Protein in sterile H₂O to a concentration of 100-200 µg/mL.
- Then, incubate it at room temperature for at least 20 mins to ensure sufficient dissolution.

Do Not Vortex! Vigorous shaking may impair the biological activity of the protein.

- Store at 2°C to 8°C for up to 1 week.

3. Extended Storage

- After the initial reconstitution, further dilute the reconstituted protein in a buffer containing a carrier protein or stabilizer (e.g., 0.1% BSA, 10% FBS, 5% HSA, or 5% trehalose solution). The final concentration is not less than 10 µg/ml.
- Prepare aliquots (≥20 µl).
- Store at -20°C or -80°C for 3 to 6 months.

Avoid repeated freeze-thaw cycles.

产品储存/运输

Lyophilized Protein: Store at -20°C for 1 year.

储存


Initial Reconstitution: Store at 2°C to 8°C for up to 1 week.

Extended Storage: Store at -20°C or -80°C for 3 to 6 months with a carrier protein or stabilizer.

运输

Blue Ice

实验数据展示



SDS-PAGE analysis of recombinant human IL-28A

For Research Use or Further Manufacturing Only

Chamot Biotechnology(Shanghai) Co., Ltd. www.chamot-bio.com

Tel: 021-51880030 Mail: info@chamot-bio.com QQ: 864920491