

# CHAMOT

**重组小鼠纤毛神经营养因子 /Recombinant CNTF,  
Mouse, Animal-Free**

CM082-5MP

CM082-20MP

CM082-100MP

CM082-500MP

CM082-1000MP



CHAMOT

乔默®生物

*Specialize In Cytokines*

# CONTENT

1 产品简介

2 产品使用

3 产品储存/运输

4 实验数据展示

## 重组小鼠纤毛神经营养因子 /Recombinant CNTF, Mouse, Animal-Free

|      |           |            |             |             |              |
|------|-----------|------------|-------------|-------------|--------------|
| 产品编号 | CM082-5MP | CM082-20MP | CM082-100MP | CM082-500MP | CM082-1000MP |
| 规格   | 5 µg      | 20 µg      | 100 µg      | 500 µg      | 1 mg         |

### 产品简介

#### 背景描述

Ciliary Neurotrophic Factor (CNTF) is the member of IL-6 cytokine family and mainly expressed in the nervous system. Mouse CNTF shares 84% sequence homology with human CNTF. CNTF is 22.9 kDa neurotrophic factor containing 110 residues, which shows multiple effects in vertebrate retinogenesis. Besides, CNTF acts as a promoter that not only accelerates adult neurogenesis but also increases the survival of neuron after injury.

#### 别称

ciliary neurotrophic factor, A1429687

#### 蛋白编码

P51642

#### 分子量

The protein has a calculated MW of 23.40 kDa. The protein migrates as 17-25 kDa under reducing condition (SDS-PAGE analysis).

#### 表达系统

Escherichia coli

#### 纯度

>98% as determined by SDS-PAGE.

#### 生物活性

Measure by its ability to induce proliferation in TF-1 cells. The ED<sub>50</sub> for this effect is <10 ng/mL. The specific activity of recombinant mouse CNTF is > 1 x 10<sup>5</sup> IU/mg.

#### 内毒素检测

<0.1 EU per 1 µg of the protein by the LAL method.

#### 蛋白序列

MAFAEQSPLTLHRRDLCSRSIWLARKIRSDLTALMESYVKHQGLNKNISLDSVDGVPVASTDR  
WSEMTEAERLQENLQAYRTFQGMLTKLLEDQRVHFTPTGDFHQAIHTLTLQVSAFAYQLEEL  
MALLEQKVPEKEADGMPVTIGDGGFLFEKKLWGLKVLQELSQWTVRSIHDLRVISSHHMGISA  
HESHYGAKQM with polyhistidine tag at the C-terminus.

#### 产品形式

The protein was lyophilized from a 0.2 µm filtered solution containing 1X PBS, pH 7.4. If you have any concerns or special requirements, please confirm with us.

#### 产品应用

Cell Culture

## 产品使用

1. Before opening, centrifuge at 3000 rpm for 5 mins.

### 2. Initial Reconstitution

- Reconstitute the Lyophilized Protein in sterile H<sub>2</sub>O to a concentration of 100-200 µg/mL.
- Then, incubate it at room temperature for at least 20 mins to ensure sufficient dissolution.

Do Not Vortex! Vigorous shaking may impair the biological activity of the protein.

- Store at 2°C to 8°C for up to 1 week.

### 3. Extended Storage

- After the initial reconstitution, further dilute the reconstituted protein in a buffer containing a carrier protein or stabilizer (e.g., 0.1% BSA, 10% FBS, 5% HSA, or 5% trehalose solution). The final concentration is not less than 10 µg/mL.

- Prepare aliquots (≥20 µl).
- Store at -20°C or -80°C for 3 to 6 months.

Avoid repeated freeze-thaw cycles.

## 产品储存/运输

**Lyophilized Protein:** Store at -20°C for 1 year.

### 储存


**Initial Reconstitution:** Store at 2°C to 8°C for up to 1 week.

**Extended Storage:** Store at -20°C or -80°C for 3 to 6 months with a carrier protein or stabilizer.

### 运输

Blue Ice

## 实验数据展示



For Research Use or Further Manufacturing Only

Chamot Biotechnology(Shanghai) Co., Ltd. [www.chamot-bio.com](http://www.chamot-bio.com)

Tel: 021-51880030 Mail: [info@chamot-bio.com](mailto:info@chamot-bio.com) QQ: 864920491