

CHAMOT

3D Cell and Organoid Culturing Gel Kit

CM001-2.5CO

CM001-5CO

CM001-10CO



CHAMOT

乔默®生物

Specialize In Cytokines



CONTENT

1 Product Specifications

2 Storage/Shipping

3D Cell and Organoid Culturing Gel Kit

Catalog#	CM001-2.5CO	CM001-5CO	CM001-10CO
Size	120 assays	240 assays	480 assays

Product Specifications

Background	The 3D cell culture hydrogel system simulates the microenvironment of human tissue using synthetic biological materials. It is stable, adjustable, easy to operate, and has high biological compatibility, which can be applied in 3D cell culture, tumor formation experiment, tumor angiogenesis and invasion study.
------------	--

Storage/Shipping

Stability & Storage	2-8°C for 24 months in lyophilized state from date of receipt.
Shipping	Blue Ice

For Research Use or Further Manufacturing Only

Chamot Biotechnology(Shanghai) Co., Ltd. www.chamot-bio.com

Tel: 021-51880030 Mail: info@chamot-bio.com QQ: 864920491