

CHAMOT

Tumor Gel Kit

CM002-2.5TX

CM002-5TX


CM002-10TX



CHAMOT

乔默®生物

Specialize In Cytokines



CONTENT

1 Product Specifications

2 Storage/Shipping

Tumor Gel Kit

Catalog#	CM002-2.5TX	CM002-5TX	CM002-10TX
Size	25 assays	50 assays	100 assays

Product Specifications

Background

Tumor Gel is made from synthetic biological materials, which can provide extracellular matrix microenvironment without the risk of zoonotic pathogens or toxins. It has high activity, stability, biocompatibility, safety, and can be operated at room temperature. It is also easy to inject and allows self-adjustment of gel density. This product can be used in tumor growth experiments with various mouse cancer models.

Storage/Shipping

Stability & Storage

2-8°C for 24 months in lyophilized state from date of receipt.

Shipping

Blue Ice

For Research Use or Further Manufacturing Only

Chamot Biotechnology(Shanghai) Co., Ltd. www.chamot-bio.com

Tel: 021-51880030 Mail: info@chamot-bio.com QQ: 864920491