
# **CHAMOT**

# **Tumor Gel Kit**

CM002-2.5TX

CM002-5TX

CM002-10TX





- 1 Product Specifications
- 2 Storage/Shipping



# **Product Datasheet**

#### **Tumor Gel Kit**

Catalog#	CM002-2.5TX	CM002-5TX	CM002-10TX
Size	25 assays	50 assays	100 assays

### **Product Specifications**

#### **Background**

Tumor Gel is made from synthetic biological materials, which can provide extracellular matrix microenvironment without the risk of zoonotic pathogens or toxins. It has high activity, stability, biocompatibility, safety, and can be operated at room temperature. It is also easy to inject and allows self-adjustment of gel density. This product can be used in tumor growth experiments with various mouse cancer models.

## Storage/Shipping

**Shipping** Blue Ice

For Research Use or Further Manufacturing Only

Chamot Biotechnology(Shanghai) Co., Ltd. www.chamot-bio.com

Tel: 021-51880030 Mail: info@chamot-bio.com QQ: 864920491