

CHAMOT

AccuBlot™ anti-Mouse IgG-HRP

CM009-0.1A2H

CM009-0.5A2H

CM009-1A2H

CM009-5A2H

CM009-10A2H



CHAMOT

乔默®生物

Specialize In Cytokines



CONTENT

1 Product Specifications

2 Product Note

3 Storage/Shipping

AccuBlot™ anti-Mouse IgG-HRP

Catalog#	CM009-0.1A2H	CM009-0.5A2H	CM009-1A2H	CM009-5A2H	CM009-10A2H
Size	100 µL	500 µL	1 mL	5 mL	10 mL

Product Specifications

Description AccuBlot™ Anti-Mouse IgG Antibody (HRP) is conjugated with high-activity HRP. In coimmunoprecipitation (co-IP) assay, denatured IgG from IP samples may show interference and result in immunoblotting at ~55 kDa (heavy chain) and ~25 kDa (light chain) by using traditional secondary antibodies. The products of AccuBlot™ series are suitable for native IgG detection.

Form Liquid

Reactivity Mouse

Clonality Monoclonal

Immunogen Native mouse IgG

Concentration Lot-dependent, please refer to product label.

Conjugation HRP

Storage Buffer Phosphate Buffered Saline containing 1% BSA, 0.03% ProClin 300 and 50% glycerol, pH 7.4.

Isotype IgG2b

Product Note

Recommended dilution factor:

ELISA: 1:2500-10000

WB:1:1000-3000

Note: Working dilution for specific application should be determined by the investigator to obtain the best conditions.

Storage/Shipping

Stability & Storage

This product is stable after storage at: 2-8°C for 2 weeks under sterile conditions from date of receipt. -20°C or -80°C for 12 months under sterile conditions from date of receipt. Avoid repeated freeze/thaw cycles and protect from light. Suggestion: Divide antibody into several vials. Keep only vials for usage at 2-8°C.

Shipping

The product is shipped with polar packs. Upon receipt, store it immediately at -20°C or lower for long term storage.

For Research Use or Further Manufacturing Only

Chamot Biotechnology(Shanghai) Co., Ltd. www.chamot-bio.com

Tel: 021-51880030 Mail: info@chamot-bio.com QQ: 864920491