

# CHAMOT

## Goat anti-Human IgG (Fc Fragment)-HRP

CM008-0.1A2H

CM008-0.5A2H

CM008-1A2H

CM008-5A2H

CM008-10A2H



CHAMOT

乔默®生物

*Specialize In Cytokines*



# CONTENT

1 Product Specifications

2 Product Note

3 Storage/Shipping

## Goat anti-Human IgG (Fc Fragment)-HRP

| Catalog# | CM008-0.1A2H | CM008-0.5A2H | CM008-1A2H | CM008-5A2H | CM008-10A2H |
|----------|--------------|--------------|------------|------------|-------------|
| Size     | 100 µL       | 500 µL       | 1 mL       | 5 mL       | 10 mL       |

### Product Specifications

|                       |                                                                                                                                                                                                   |
|-----------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Description</b>    | Goat anti-Human IgG, Fc Fragment is purified from goat antiserum by immunoaffinity chromatography and conjugated with high-activity HRP. Minimal cross-reactivity with other isotypes or species. |
| <b>Form</b>           | Liquid                                                                                                                                                                                            |
| <b>Reactivity</b>     | Human                                                                                                                                                                                             |
| <b>Clonality</b>      | Polyclonal                                                                                                                                                                                        |
| <b>Immunogen</b>      | IgG, Fc Fragment                                                                                                                                                                                  |
| <b>Concentration</b>  | Lot-dependent, please refer to product label.                                                                                                                                                     |
| <b>Conjugation</b>    | HRP                                                                                                                                                                                               |
| <b>Storage Buffer</b> | Phosphate Buffered Saline containing 1% BSA, 0.03% ProClin 300, pH 7.4.                                                                                                                           |
| <b>Isotype</b>        | IgG, Fc Fragment                                                                                                                                                                                  |

### Product Note

Recommended dilution factor:

ELISA: 1:5000-25000

IHC: 1:500-2500

WB:1:5000-25000

Note: Working dilution for specific application should be determined by the investigator to obtain the best conditions.

## Storage/Shipping

### Stability & Storage

This product is stable after storage at: 2-8°C for 2 weeks under sterile conditions from date of receipt. -20°C or -80°C for 12 months under sterile conditions from date of receipt. Avoid repeated freeze/thaw cycles and protect from light. Suggestion: Divide antibody into several vials. Keep only vials for usage at 2-8°C.

### Shipping

The product is shipped with polar packs. Upon receipt, store it immediately at -20°C or lower for long term storage.

For Research Use or Further Manufacturing Only

Chamot Biotechnology(Shanghai) Co., Ltd. [www.chamot-bio.com](http://www.chamot-bio.com)

Tel: 021-51880030    Mail: [info@chamot-bio.com](mailto:info@chamot-bio.com)    QQ: 864920491