

# CHAMOT

## Goat anti-Human IgG, Fc Fragment

CM007-0.1A2

CM007-0.5A2

CM007-1A2

CM007-5A2

CM007-10A2



CHAMOT

乔默®生物

*Specialize In Cytokines*



# CONTENT

1 Product Specifications

2 Storage/Shipping

## Goat anti-Human IgG, Fc Fragment

|          |             |             |           |           |            |
|----------|-------------|-------------|-----------|-----------|------------|
| Catalog# | CM007-0.1A2 | CM007-0.5A2 | CM007-1A2 | CM007-5A2 | CM007-10A2 |
| Size     | 100 µg      | 500 µg      | 1 mg      | 5 mg      | 10 mg      |

### Product Specifications

|                       |                                                                                                                                                            |
|-----------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Description</b>    | Goat anti-Human IgG, Fc Fragment is purified from goat antiserum by immunoaffinity chromatography. Minimal cross-reactivity with other isotype or species. |
| <b>Form</b>           | Liquid                                                                                                                                                     |
| <b>Reactivity</b>     | Human                                                                                                                                                      |
| <b>Clonality</b>      | Polyclonal                                                                                                                                                 |
| <b>Immunogen</b>      | Gamma immunoglobulin Fc region                                                                                                                             |
| <b>Concentration</b>  | 1 mg/mL                                                                                                                                                    |
| <b>Conjugation</b>    | Unconjugated                                                                                                                                               |
| <b>Storage Buffer</b> | Phosphate Buffered Saline containing 0.05% sodium azide, pH 7.4.                                                                                           |
| <b>Isotype</b>        | IgG, Fc Fragment                                                                                                                                           |

### Storage/Shipping

|                                |                                                                                                                                                                                                                                                                                                              |
|--------------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Stability &amp; Storage</b> | This product is stable after storage at: 2-8°C for 2 weeks under sterile conditions from date of receipt. -20°C or -80°C for 12 months under sterile conditions from date of receipt. Avoid repeated freeze/thaw cycles. Suggestion: Divide antibody into several vials. Keep only vials for usage at 2-8°C. |
| <b>Shipping</b>                | The product is shipped with polar packs. Upon receipt, store it immediately at -20°C or lower for long term storage.                                                                                                                                                                                         |

For Research Use or Further Manufacturing Only

Chamot Biotechnology(Shanghai) Co., Ltd. [www.chamot-bio.com](http://www.chamot-bio.com)

Tel: 021-51880030    Mail: [info@chamot-bio.com](mailto:info@chamot-bio.com)    QQ: 864920491