

# CHAMOT

## Goat anti-Mouse IgG (H+L)

CM002-0.1A2

CM002-0.5A2

CM002-1A2

CM002-5A2

CM002-10A2



CHAMOT

乔默®生物

*Specialize In Cytokines*



# CONTENT

1 Product Specifications

2 Storage/Shipping

## Goat anti-Mouse IgG (H+L)

Catalog#	CM002-0.1A2	CM002-0.5A2	CM002-1A2	CM002-5A2	CM002-10A2
Size	100 µg	500 µg	1 mg	5 mg	10 mg

### Product Specifications

<b>Description</b>	Goat Anti-Mouse IgG (H+L) Antibody is purified from goat, specifically react with mouse IgG, heavy and light chain.
<b>Form</b>	Liquid
<b>Reactivity</b>	Mouse (Minor cross-reactivity to rat)
<b>Clonality</b>	Polyclonal
<b>Immunogen</b>	Purified mouse IgG
<b>Concentration</b>	5 mg/mL
<b>Conjugation</b>	Unconjugated
<b>Storage Buffer</b>	Phosphate Buffered Saline containing 0.05% sodium azide, pH 7.4.
<b>Isotype</b>	IgG

### Storage/Shipping

<b>Stability &amp; Storage</b>	This product is stable after storage at: 2-8°C for 2 weeks under sterile conditions from date of receipt. -20°C or -80°C for 12 months under sterile conditions from date of receipt. Avoid repeated freeze/thaw cycles. Suggestion: Divide antibody into several vials. Keep only vials for usage at 2-8°C.
<b>Shipping</b>	The product is shipped with polar packs. Upon receipt, store it immediately at -20°C or lower for long term storage.

For Research Use or Further Manufacturing Only

Chamot Biotechnology(Shanghai) Co., Ltd. [www.chamot-bio.com](http://www.chamot-bio.com)

Tel: 021-51880030    Mail: [info@chamot-bio.com](mailto:info@chamot-bio.com)    QQ: 864920491