

CHAMOT

CM011-0.1A2F

CM011-0.5A2F

CM011-1A2F

CM011-5A2F

CM011-10A2F



CHAMOT

乔默®生物

Specialize In Cytokines



CONTENT

1 Product Specifications

2 Product Note

3 Storage/Shipping

Goat anti-Mouse IgG-iFluor 405

Catalog#	CM011-0.1A2F	CM011-0.5A2F	CM011-1A2F	CM011-5A2F	CM011-10A2F
Size	100 µL	500 µL	1 mL	5 mL	10 mL

Product Specifications

Description iFluor 405 is a blue-fluorescent dye that can be excited by the 405 nm laser. LeadgeneGoat Anti-Mouse IgG (H+L) Antibody (iFluor 405) is conjugated with high-activity iFluor 405.

Form Liquid

Reactivity Mouse (Minor cross-reactivity to rat)

Clonality Polyclonal

Immunogen Purified mouse IgG

Concentration 0.5 mg/mL

Conjugation iFluor 405

Storage Buffer Phosphate Buffered Saline containing 1% BSA, 0.03% ProClin 300, pH 7.4.

Isotype IgG

Product Note

Recommended dilution factor:

IF: 1:200-1000

FACS:1:200-1000

Note: Working dilution for specific application should be determined by the investigator to obtain the best conditions.

Storage/Shipping

Stability & Storage

This product is stable after storage at: 2-8°C for 6 months under sterile conditions from date of receipt. Do not freeze and protect from light.

Shipping

The product is shipped with polar packs. Upon receipt, store it immediately at 2-8°C for long term storage.

For Research Use or Further Manufacturing Only

Chamot Biotechnology(Shanghai) Co., Ltd. www.chamot-bio.com

Tel: 021-51880030 Mail: info@chamot-bio.com QQ: 864920491