

# CHAMOT

CM011-0.1A2F

CM011-0.5A2F

CM011-1A2F

CM011-5A2F


CM011-10A2F



CHAMOT

乔默®生物

Specialize In Cytokines



# CONTENT

1 Product Specifications

2 Product Note

3 Storage/Shipping

## Goat anti-Mouse IgG-iFluor 405

Catalog#	CM011-0.1A2F	CM011-0.5A2F	CM011-1A2F	CM011-5A2F	CM011-10A2F
Size	100 µL	500 µL	1 mL	5 mL	10 mL

### Product Specifications

**Description** iFluor 405 is a blue-fluorescent dye that can be excited by the 405 nm laser. LeadgeneGoat Anti-Mouse IgG (H+L) Antibody (iFluor 405) is conjugated with high-activity iFluor 405.

**Form** Liquid

**Reactivity** Mouse (Minor cross-reactivity to rat)

**Clonality** Polyclonal

**Immunogen** Purified mouse IgG

**Concentration** 0.5 mg/mL

**Conjugation** iFluor 405

**Storage Buffer** Phosphate Buffered Saline containing 1% BSA, 0.03% ProClin 300, pH 7.4.

**Isotype** IgG

### Product Note

Recommended dilution factor:

IF: 1:200-1000

FACS:1:200-1000

Note: Working dilution for specific application should be determined by the investigator to obtain the best conditions.

## Storage/Shipping

### **Stability & Storage**

This product is stable after storage at: 2-8°C for 6 months under sterile conditions from date of receipt. Do not freeze and protect from light.

### **Shipping**

The product is shipped with polar packs. Upon receipt, store it immediately at 2-8°C for long term storage.

For Research Use or Further Manufacturing Only

Chamot Biotechnology(Shanghai) Co., Ltd. [www.chamot-bio.com](http://www.chamot-bio.com)

Tel: 021-51880030    Mail: [info@chamot-bio.com](mailto:info@chamot-bio.com)    QQ: 864920491