

CHAMOT

CM013-0.1A2F

CM013-0.5A2F

CM013-1A2F

CM013-5A2F

CM013-10A2F



CHAMOT

乔默®生物

Specialize In Cytokines



CONTENT

1 Product Specifications

2 Product Note

3 Storage/Shipping

Goat anti-Mouse IgG-iFluor 594

Catalog#	CM013-0.1A2F	CM013-0.5A2F	CM013-1A2F	CM013-5A2F	CM013-10A2F
Size	100 µL	500 µL	1 mL	5 mL	10 mL

Product Specifications

Description	iFluor 594 dye is a bright, red-fluorescent dye that can be excited by the 561 nm or 594 nm laser. Leadgene Goat Anti-Mouse IgG (H+L) Antibody (iFluor 594) is conjugated with high-activity iFluor 594.
Form	Liquid
Reactivity	Mouse (Minor cross-reactivity to rat)
Clonality	Polyclonal
Immunogen	Purified mouse IgG
Concentration	0.5 mg/mL
Conjugation	iFluor 594
Storage Buffer	Phosphate Buffered Saline containing 1% BSA, 0.03% ProClin 300, pH 7.4.
Isotype	IgG

Product Note

Recommended dilution factor:

IF: 1:200-1000

FACS:1:200-1000

Note: Working dilution for specific application should be determined by the investigator to obtain the best conditions.

Storage/Shipping

Stability & Storage

This product is stable after storage at: 2-8°C for 6 months under sterile conditions from date of receipt. Do not freeze and protect from light.

Shipping

The product is shipped with polar packs. Upon receipt, store it immediately at 2-8°C for long term storage.

For Research Use or Further Manufacturing Only

Chamot Biotechnology(Shanghai) Co., Ltd. www.chamot-bio.com

Tel: 021-51880030 Mail: info@chamot-bio.com QQ: 864920491