

# CHAMOT

CM003-0.1A2F

CM003-0.5A2F

CM003-1A2F

CM003-5A2F

CM003-10A2F



CHAMOT

乔默®生物

*Specialize In Cytokines*



# CONTENT

1 Product Specifications

2 Product Note

3 Storage/Shipping

## Rabbit anti-Mouse IgG-iFluor 594

|          |              |              |            |            |             |
|----------|--------------|--------------|------------|------------|-------------|
| Catalog# | CM003-0.1A2F | CM003-0.5A2F | CM003-1A2F | CM003-5A2F | CM003-10A2F |
| Size     | 100 µL       | 500 µL       | 1 mL       | 5 mL       | 10 mL       |

### Product Specifications

|                       |                                                                                                                                                                                                            |
|-----------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Description</b>    | iFluor 594 dye is a bright, red-fluorescent dye that can be excited by the 561 nm or 594 nm laser. Leadgene Rabbit Anti-Mouse IgG (H+L) Antibody (iFluor 594) is conjugated with high-activity iFluor 594. |
| <b>Form</b>           | Liquid                                                                                                                                                                                                     |
| <b>Reactivity</b>     | Mouse                                                                                                                                                                                                      |
| <b>Clonality</b>      | Polyclonal                                                                                                                                                                                                 |
| <b>Immunogen</b>      | Purified mouse IgG                                                                                                                                                                                         |
| <b>Concentration</b>  | 0.5 mg/mL                                                                                                                                                                                                  |
| <b>Conjugation</b>    | iFluor 594                                                                                                                                                                                                 |
| <b>Storage Buffer</b> | Phosphate Buffered Saline containing 1% BSA, 0.03% ProClin 300, pH 7.4.                                                                                                                                    |
| <b>Isotype</b>        | IgG                                                                                                                                                                                                        |

### Product Note

Recommended dilution factor:

IF: 1:200-1000

FACS:1:200-1000

Note: Working dilution for specific application should be determined by the investigator to obtain the best conditions.

## Storage/Shipping

**Stability & Storage**

This product is stable after storage at: 2-8°C for 6 months under sterile conditions from date of receipt. Do not freeze and protect from light.

**Shipping**

The product is shipped with polar packs. Upon receipt, store it immediately at 2-8°C for long term storage.

For Research Use or Further Manufacturing Only

Chamot Biotechnology(Shanghai) Co., Ltd. [www.chamot-bio.com](http://www.chamot-bio.com)

Tel: 021-51880030    Mail: [info@chamot-bio.com](mailto:info@chamot-bio.com)    QQ: 864920491