

CHAMOT

Mouse anti-DYKDDDDK tag mAb

CM006-0.1A1

CM006-0.5A1

CM006-1A1

CM006-5A1

CM006-10A1



CHAMOT

乔默®生物

Specialize In Cytokines



CONTENT

1 Product Specifications

2 Product Note

3 Storage/Shipping

Mouse anti-DYKDDDDK tag mAb

Catalog#	CM006-0.1A1	CM006-0.5A1	CM006-1A1	CM006-5A1	CM006-10A1
Size	100 µg	500 µg	1 mg	5 mg	10 mg

Product Specifications

Description	Mouse anti-DYKDDDDK tag mAb is specific to DYKDDDDK (FLAG) tag, which can be applied to several immunoassays, including ELISA, WB, IP, IFA, and IHC. This antibody has been identified in several sources (bacteria, yeast, and mammalian cells) of FLAG tag fusion protein.
Form	Liquid
Reactivity	Recognizes tag sequence of DYKDDDDK on either the N- or C-terminus of recombinant fusion protein.
Clonality	Monoclonal, clone 108-21
Immunogen	Pepetide: DYKDDDDK (FLAG)
Concentration	0.5 mg/mL
Conjugation	Unconjugated
Storage Buffer	Phosphate Buffered Saline containing 50% glycerol and 0.02% sodium azide.
Isotype	IgG2b

Product Note

Recommended dilution factor:

ELISA: 1:5000-20000

WB:1:1000-10000

IFA: 1:500-2000

IHC: 1:500-2000

IP: 1:100

Note: Working dilution for specific application should be determined by the investigator to obtain the best conditions.

Storage/Shipping

Stability & Storage

This product is stable after storage at: 2-8°C for 2 weeks under sterile conditions from date of receipt. -20°C or -80°C for 12 months under sterile conditions from date of receipt. Avoid repeated freeze/thaw cycles. Suggestion: Divide antibody into several vials. Keep only vials for usage at 2-8°C.

Shipping

The product is shipped with polar packs. Upon receipt, store it immediately at -20°C or lower for long term storage.

For Research Use or Further Manufacturing Only

Chamot Biotechnology(Shanghai) Co., Ltd. www.chamot-bio.com

Tel: 021-51880030 Mail: info@chamot-bio.com QQ: 864920491