

CHAMOT

Mouse anti-DYKDDDDK tag mAb

CM006-A1S

CM006-0.1A1

CM006-0.5A1

CM006-1A1

CM006-5A1

CM006-10A1



CHAMOT

乔默®生物

Specialize In Cytokines



CONTENT

1 产品简介

2 产品使用

3 产品储存/运输

Mouse anti-DYKDDDDK tag mAb

产品编号	CM006-A1S	CM006-0.1A1	CM006-0.5A1	CM006-1A1	CM006-5A1	CM006-10A1
规格	10 µg	100 µg	500 µg	1 mg	5 mg	10 mg

产品简介

产品描述

Mouse anti-DYKDDDDK tag mAb is specific to DYKDDDDK (FLAG) tag, which can be applied to several immunoassays, including ELISA, WB, IP, IFA, and IHC. This antibody has been identified in several sources (bacteria, yeast, and mammalian cells) of FLAG tag fusion protein.

产品形式

Liquid

产品应用

ELISA,WB,IFA,IHC,IP

反应性

Recognizes tag sequence of DYKDDDDK on either the N- or C-terminus of recombinant fusion protein.

克隆性

Monoclonal, clone 108-21

免疫原

Pepetide: DYKDDDDK (FLAG)

浓度

0.5 mg/mL

标记

Unconjugated

产品组成

Phosphate Buffered Saline containing 50% glycerol and 0.02% sodium azide.

亚型

IgG2b

产品使用

Recommended dilution factor:

ELISA: 1:5000-20000

WB:1:1000-10000

IFA: 1:500-2000

IHC: 1:500-2000

IP: 1:100

Note: Working dilution for specific application should be determined by the investigator to obtain the best conditions.

产品储存/运输

储存

This product is stable after storage at: 2-8°C for 2 weeks under sterile conditions from date of receipt. -20°C or -80°C for 12 months under sterile conditions from date of receipt. Avoid repeated freeze/thaw cycles. Suggestion: Divide antibody into several vials. Keep only vials for usage at 2-8°C.

运输

The product is shipped with polar packs. Upon receipt, store it immediately at -20°C or lower for long term storage.

For Research Use or Further Manufacturing Only

Chamot Biotechnology(Shanghai) Co., Ltd. www.chamot-bio.com

Tel: 021-51880030 Mail: info@chamot-bio.com QQ: 864920491