

CHAMOT

Mouse anti-GST tag mAb

CM005-0.1A1

CM005-0.5A1

CM005-1A1

CM005-5A1

CM005-10A1



CHAMOT

乔默®生物

Specialize In Cytokines



CONTENT

1 Product Specifications

2 Product Note

3 Storage/Shipping

Mouse anti-GST tag mAb

Catalog#	CM005-0.1A1	CM005-0.5A1	CM005-1A1	CM005-5A1	CM005-10A1
Size	100 µg	500 µg	1 mg	5 mg	10 mg

Product Specifications

Description	Mouse anti-GST tag mAb is specific for native as well as denatured-reduced forms of GST (glutathione S-transferase) tag, which can be applied to several immunoassays, including ELISA, WB, IP, IFA, and FACS.
Form	Liquid
Reactivity	Recognizes glutathione S-transferase tag on either the N- or C-terminus of recombinant fusion protein.
Clonality	Monoclonal, clone 89-3
Immunogen	Recombinant glutathione S-transferase
Concentration	1 mg/mL
Conjugation	Unconjugated
Storage Buffer	Phosphate Buffered Saline containing 0.03% ProClin 300, pH 7.4.
Isotype	IgG2a

Product Note

Recommended dilution factor:

ELISA: 1:5000-20000

WB: 1:1000-5000

IP: 1:500-2000

IFA: 1:500-1000

FACS: 1:500-1000

Note: Working dilution for specific application should be determined by the investigator to obtain the best conditions.

Storage/Shipping

Stability & Storage

This product is stable after storage at: 2-8°C for 2 weeks under sterile conditions from date of receipt. -20°C or -80°C for 12 months under sterile conditions from date of receipt. Avoid repeated freeze/thaw cycles. Suggestion: Divide antibody into several vials. Keep only vials for usage at 2-8°C.

Shipping

The product is shipped with polar packs. Upon receipt, store it immediately at -20°C or lower for long term storage.

For Research Use or Further Manufacturing Only

Chamot Biotechnology(Shanghai) Co., Ltd. www.chamot-bio.com

Tel: 021-51880030 Mail: info@chamot-bio.com QQ: 864920491