

CHAMOT

Mouse anti-HA tag mAb

CM001-0.1A1

CM001-0.5A1

CM001-1A1

CM001-5A1

CM001-10A1



CHAMOT

乔默®生物

Specialize In Cytokines



CONTENT

1 Product Specifications

2 Product Note

3 Storage/Shipping

Mouse anti-HA tag mAb

Catalog#	CM001-0.1A1	CM001-0.5A1	CM001-1A1	CM001-5A1	CM001-10A1
Size	100 µg	500 µg	1 mg	5 mg	10 mg

Product Specifications

Description	Mouse anti-HA tag mAb is specific to YPYDVPDYA (HA) tag which can be applied to several immunoassays including ELISA, WB, IP, IFA and FACS. This antibody has been identified in several sources (bacteria, yeast and mammalian cells) of HA tag fusion protein.
Form	Liquid
Reactivity	Recognizes HA tag on either the N- or C-terminus of recombinant fusion protein.
Clonality	Monoclonal, clone 98-106
Immunogen	Peptide: YPYDVPDYA (HA)
Concentration	1 mg/mL
Conjugation	Unconjugated
Storage Buffer	Phosphate Buffered Saline containing 50% glycerol and 0.02% sodium azide.

Product Note

Recommended dilution factor:

ELISA: 1:5000-20000

WB: 1:1000-5000

IP: 1:500-2000

IFA: 1:500-1000

FACS: 1:500-1000

Note: Working dilution for specific application should be determined by the investigator to obtain the best conditions.

Storage/Shipping

Stability & Storage

This product is stable after storage at: 2-8°C for 2 weeks under sterile conditions from date of receipt. -20°C or -80°C for 12 months under sterile conditions from date of receipt. Avoid repeated freeze/thaw cycles. Suggestion: Divide antibody into several vials. Keep only vials for usage at 2-8°C.

Shipping

The product is shipped with polar packs. Upon receipt, store it immediately at -20°C or lower for long term storage.

For Research Use or Further Manufacturing Only

Chamot Biotechnology(Shanghai) Co., Ltd. www.chamot-bio.com

Tel: 021-51880030 Mail: info@chamot-bio.com QQ: 864920491