

CHAMOT

Mouse anti-His tag mAb

CM002-0.1A1

CM002-0.5A1

CM002-1A1

CM002-5A1

CM002-10A1



CHAMOT

乔默®生物

Specialize In Cytokines



CONTENT

1 Product Specifications

2 Product Note

3 Storage/Shipping

Mouse anti-His tag mAb

Catalog#	CM002-0.1A1	CM002-0.5A1	CM002-1A1	CM002-5A1	CM002-10A1
Size	100 µg	500 µg	1 mg	5 mg	10 mg

Product Specifications

Description Mouse anti-His tag mAb is specific to 6X His tag, which can be applied to several immunoassays, including ELISA, WB, IP, IFA, and FACS. This antibody has been identified in several sources (bacteria, yeast, and mammalian cells) of His tag fusion protein.

Form Liquid

Reactivity Recognizes 6X His tag on either the N- or C-terminus of recombinant fusion protein.

Clonality Monoclonal, clone 14-4

Immunogen 6X His tag fusion protein

Concentration 1 mg/mL

Conjugation Unconjugated

Storage Buffer Phosphate Buffered Saline containing 0.03% ProClin 300, pH 7.4.

Isotype IgG1

Product Note

Recommended dilution factor:

ELISA: 1:5000-20000

WB: 1:1000-5000

IP: 1:500-2000

IFA: 1:500-1000

FACS: 1:500-1000

Note: Working dilution for specific application should be determined by the investigator to obtain the best conditions.

Storage/Shipping

Stability & Storage

This product is stable after storage at: 2-8°C for 2 weeks under sterile conditions from date of receipt. -20°C or -80°C for 12 months under sterile conditions from date of receipt. Avoid repeated freeze/thaw cycles. Suggestion: Divide antibody into several vials. Keep only vials for usage at 2-8°C.

Shipping

The product is shipped with polar packs. Upon receipt, store it immediately at -20°C or lower for long term storage.

For Research Use or Further Manufacturing Only

Chamot Biotechnology(Shanghai) Co., Ltd. www.chamot-bio.com

Tel: 021-51880030 Mail: info@chamot-bio.com QQ: 864920491