

CHAMOT

Mouse anti-mCherry mAb

CM008-0.1A1

CM008-0.5A1

CM008-1A1

CM008-5A1

CM008-10A1



CHAMOT

乔默®生物

Specialize In Cytokines



CONTENT

1 Product Specifications

2 Product Note

3 Storage/Shipping

Mouse anti-mCherry mAb

Catalog#	CM008-0.1A1	CM008-0.5A1	CM008-1A1	CM008-5A1	CM008-10A1
Size	100 µg	500 µg	1 mg	5 mg	10 mg

Product Specifications

Description	The mCherry protein is a basic red fluorescent protein derived from <i>Discosoma</i> sp. Leadgene mouse anti-mCherry mAb can recognize native and denature mCherry protein but not GFP.
Form	Liquid
Reactivity	Recognizes mCherry but not GFP.
Clonality	Monoclonal, clone 45-5
Immunogen	Recombinant mCherry protein
Concentration	1 mg/mL
Conjugation	Unconjugated
Storage Buffer	Phosphate Buffered Saline containing 0.03% ProClin 300, pH 7.4.
Isotype	IgG2a

Product Note

Recommended dilution factor:

ELISA: 1:5000-20000

WB: 1:1000-5000

IP: 1:500-2000

IFA: 1:500-1000

FACS: 1:500-1000

Note: Working dilution for specific application should be determined by the investigator to obtain the best conditions.

Storage/Shipping

Stability & Storage

This product is stable after storage at: 2-8°C for 2 weeks under sterile conditions from date of receipt. -20°C or -80°C for 12 months under sterile conditions from date of receipt. Avoid repeated freeze/thaw cycles. Suggestion: Divide antibody into several vials. Keep only vials for usage at 2-8°C.

Shipping

The product is shipped with polar packs. Upon receipt, store it immediately at -20°C or lower for long term storage.

For Research Use or Further Manufacturing Only

Chamot Biotechnology(Shanghai) Co., Ltd. www.chamot-bio.com

Tel: 021-51880030 Mail: info@chamot-bio.com QQ: 864920491