

CHAMOT

Mouse anti-V5 tag mAb

CM004-A1S

CM004-0.1A1

CM004-0.5A1

CM004-1A1

CM004-5A1

CM004-10A1



CHAMOT

乔默®生物

Specialize In Cytokines



CONTENT

1 产品简介

2 产品使用

3 产品储存/运输

Mouse anti-V5 tag mAb

产品编号	CM004-A1S	CM004-0.1A1	CM004-0.5A1	CM004-1A1	CM004-5A1	CM004-10A1
规格	10 µg	100 µg	500 µg	1 mg	5 mg	10 mg

产品简介

产品描述

Mouse monoclonal anti-V5 antibody is an antibody specific to GKIPNPLLGLDST (V5) tags which can be applied to several immunoassays including ELISA, WB, IHC, and Flow cytometry. This antibody has been identified in several sources (bacteria, yeast, and mammalian cells) of V5-fusion protein.

产品形式

Liquid

产品应用

ELISA, WB, IP, IFA, FACS

反应性

Recognizes tag sequence of GKIPNPLLGLDST on either the N- or C-terminus of recombinant fusion protein.

克隆性

Monoclonal, clone 7-1-8

免疫原

Peptide: GKIPNPLLGLDST (V5)

浓度

1 mg/mL

标记

Unconjugated

产品组成

Phosphate Buffered Saline containing 0.03% ProClin 300, pH 7.4.

亚型

IgG2b

产品使用

Recommended dilution factor:

ELISA: 1:5000-20000

WB: 1:1000-5000

IP: 1:500-2000

IFA: 1:500-1000
FACS: 1:500-1000

Note: Working dilution for specific application should be determined by the investigator to obtain the best conditions.

产品储存/运输

储存

This product is stable after storage at: 2-8°C for 2 weeks under sterile conditions from date of receipt. -20°C or -80°C for 12 months under sterile conditions from date of receipt. Avoid repeated freeze/thaw cycles. Suggestion: Divide antibody into several vials. Keep only vials for usage at 2-8°C.

运输

The product is shipped with polar packs. Upon receipt, store it immediately at -20°C or lower for long term storage.

For Research Use or Further Manufacturing Only

Chamot Biotechnology(Shanghai) Co., Ltd. www.chamot-bio.com

Tel: 021-51880030 Mail: info@chamot-bio.com QQ: 864920491