

# CHAMOT

CM007-0.1A2H

CM007-0.5A2H

CM007-1A2H

CM007-5A2H

CM007-10A2H



CHAMOT

乔默®生物

*Specialize In Cytokines*



# CONTENT

1 Product Specifications

2 Product Note

3 Storage/Shipping

## AccuBlot™ anti-Rabbit IgG-HRP

Catalog#	CM007-0.1A2H	CM007-0.5A2H	CM007-1A2H	CM007-5A2H	CM007-10A2H
Size	100 µL	500 µL	1 mL	5 mL	10 mL

### Product Specifications

**Description** AccuBlot™ Anti-Rabbit IgG (H+L) Antibody (HRP) is conjugated with high-activity HRP. In coimmunoprecipitation (co-IP) assay, denatured IgG from IP samples may show interference and result in immunoblotting at ~55 kDa (heavy chain) and ~25 kDa (light chain) by using traditional secondary antibodies. The products of AccuBlot™ series are suitable for native IgG detection.

**Form** Liquid

**Reactivity** Rabbit

**Clonality** Monoclonal

**Immunogen** Purified rabbit IgG

**Concentration** Lot-dependent, please refer to product label.

**Conjugation** HRP

**Storage Buffer** Phosphate Buffered Saline containing 1% BSA, 0.03% ProClin 300 and 50% glycerol, pH 7.4.

**Isotype** IgG2b

### Product Note

Recommended dilution factor:

ELISA: 1:2500-10000

WB:1:1000-3000

Note: Working dilution for specific application should be determined by the investigator to obtain the best conditions.

## Storage/Shipping

### Stability & Storage

This product is stable after storage at: 2-8°C for 2 weeks under sterile conditions from date of receipt. -20°C or -80°C for 12 months under sterile conditions from date of receipt. Avoid repeated freeze/thaw cycles and protect from light. Suggestion: Divide antibody into several vials. Keep only vials for usage at 2-8°C.

### Shipping

The product is shipped with polar packs. Upon receipt, store it immediately at -20°C or lower for long term storage.

For Research Use or Further Manufacturing Only

Chamot Biotechnology(Shanghai) Co., Ltd. [www.chamot-bio.com](http://www.chamot-bio.com)

Tel: 021-51880030    Mail: [info@chamot-bio.com](mailto:info@chamot-bio.com)    QQ: 864920491