

CHAMOT

AccuBlot™ anti-Rabbit IgG-HRP

CM007-A2HS

CM007-0.1A2H

CM007-0.5A2H

CM007-1A2H

CM007-5A2H


CM007-10A2H



CHAMOT

乔默®生物

Specialize In Cytokines



CONTENT

1 产品简介

2 产品使用

3 产品储存/运输

AccuBlot™ anti-Rabbit IgG-HRP

产品编号	CM007-A2HS	CM007-0.1A2H	CM007-0.5A2H	CM007-1A2H	CM007-5A2H	CM007-10A2H
规格	10 µL	100 µL	500 µL	1 mL	5 mL	10 mL

产品简介

产品描述 AccuBlot™ Anti-Rabbit IgG (H+L) Antibody (HRP) is conjugated with high-activity HRP. In coimmunoprecipitation (co-IP) assay, denatured IgG from IP samples may show interference and result in immunoblotting at ~55 kDa (heavy chain) and ~25 kDa (light chain) by using traditional secondary antibodies. The products of AccuBlot™ series are suitable for native IgG detection.

产品形式 Liquid

产品应用 ELISA, WB

反应性 Rabbit IgG

克隆性 Monoclonal

免疫原 Purified rabbit IgG

浓度 Lot-dependent, please refer to product label.

标记 HRP

产品组成 Phosphate Buffered Saline containing 1% BSA, 0.03% ProClin 300 and 50% glycerol, pH 7.4.

产品使用

Recommended dilution factor:

ELISA: 1:2500-10000

WB: 1:1000-3000

Note: Working dilution for specific application should be determined by the investigator to obtain the best conditions.

产品储存/运输

储存

This product is stable after storage at: 2-8°C for 2 weeks under sterile conditions from date of receipt. -20°C or -80°C for 12 months under sterile conditions from date of receipt. Avoid repeated freeze/thaw cycles and protect from light. Suggestion: Divide antibody into several vials. Keep only vials for usage at 2-8°C.

运输

The product is shipped with polar packs. Upon receipt, store it immediately at -20°C or lower for long term storage.

For Research Use or Further Manufacturing Only

Chamot Biotechnology(Shanghai) Co., Ltd. www.chamot-bio.com

Tel: 021-51880030 Mail: info@chamot-bio.com QQ: 864920491