

# CHAMOT

## Goat anti-Rabbit IgG-iFluor 405

CM004-A2FS

CM004-0.1A2F

CM004-0.5A2F

CM004-1A2F

CM004-5A2F

CM004-10A2F



CHAMOT

乔默®生物

*Specialize In Cytokines*



# CONTENT

1 产品简介

2 产品使用

3 产品储存/运输

## Goat anti-Rabbit IgG-iFluor 405

|      |            |              |              |            |            |             |
|------|------------|--------------|--------------|------------|------------|-------------|
| 产品编号 | CM004-A2FS | CM004-0.1A2F | CM004-0.5A2F | CM004-1A2F | CM004-5A2F | CM004-10A2F |
| 规格   | 10 µL      | 100 µL       | 500 µL       | 1 mL       | 5 mL       | 10 mL       |

### 产品简介

**产品描述** iFluor 405 is a blue-fluorescent dye that can be excited by the 405 nm laser. Leadgene Goat Anti-Rabbit IgG (H+L) Antibody (iFluor 405) is conjugated with high-activity iFluor 405.

**产品形式** Liquid

**产品应用** IF,FACS

**反应性** Rabbit (Minor cross-reactivity to human in higher concentration)

**克隆性** Polyclonal

**免疫原** Purified rabbit IgG

**浓度** 0.5 mg/mL

**标记** iFluor 405

**产品组成** Phosphate Buffered Saline containing 1% BSA, 0.03% ProClin 300, pH 7.4.

**亚型** IgG

### 产品使用

Recommended dilution factor:

IF: 1:200-1000

FACS:1:200-1000

Note: Working dilution for specific application should be determined by the investigator to obtain the best conditions.

## 产品储存/运输

### 储存

This product is stable after storage at: 2-8°C for 6 months under sterile conditions from date of receipt. Do not freeze and protect from light.

### 运输

The product is shipped with polar packs. Upon receipt, store it immediately at 2-8°C for long term storage.

For Research Use or Further Manufacturing Only

Chamot Biotechnology(Shanghai) Co., Ltd. [www.chamot-bio.com](http://www.chamot-bio.com)

Tel: 021-51880030    Mail: [info@chamot-bio.com](mailto:info@chamot-bio.com)    QQ: 864920491