

CHAMOT

Goat anti-Rabbit IgG (H+L)

CM003-A2S

CM003-0.1A2

CM003-0.5A2

CM003-1A2

CM003-5A2

CM003-10A2



CHAMOT

乔默®生物

Specialize In Cytokines



CONTENT

1 产品简介

2 产品储存/运输

Goat anti-Rabbit IgG (H+L)

产品编号	CM003-A2S	CM003-0.1A2	CM003-0.5A2	CM003-1A2	CM003-5A2	CM003-10A2
规格	10 µg	100 µg	500 µg	1 mg	5 mg	10 mg

产品简介

产品描述	Goat Anti-Rabbit IgG (H+L) Antibody is purified from goat, specifically react with rabbit IgG, heavy and light chain.
产品形式	Liquid
反应性	Rabbit (Minor cross-reactivity to human in higher concentration)
克隆性	Polyclonal
免疫原	Purified rabbit IgG
浓度	Lot-dependent, please refer to product label.
标记	Unconjugated
产品组成	Phosphate Buffered Saline containing 0.05% sodium azide, pH 7.4.
亚型	IgG

产品储存/运输

储存	This product is stable after storage at: 2-8°C for 2 weeks under sterile conditions from date of receipt. -20°C or -80°C for 12 months under sterile conditions from date of receipt. Avoid repeated freeze/thaw cycles. Suggestion: Divide antibody into several vials. Keep only vials for usage at 2-8°C.
运输	The product is shipped with polar packs. Upon receipt, store it immediately at -20°C or lower for long term storage.

For Research Use or Further Manufacturing Only

Chamot Biotechnology(Shanghai) Co., Ltd. www.chamot-bio.com

Tel: 021-51880030 Mail: info@chamot-bio.com QQ: 864920491