

CHAMOT

Mouse anti-EGFP mAb, clone 42-2

CM007-A1S

CM007-0.1A1

CM007-0.5A1

CM007-1A1

CM007-5A1

CM007-10A1



CHAMOT

乔默®生物

Specialize In Cytokines



CONTENT

1 产品简介

2 产品使用

3 产品储存/运输

Mouse anti-EGFP mAb, clone 42-2

产品编号	CM007-A1S	CM007-0.1A1	CM007-0.5A1	CM007-1A1	CM007-5A1	CM007-10A1
规格	10 µg	100 µg	500 µg	1 mg	5 mg	10 mg

产品简介

产品描述	EGFP is a 27 kDa fluorescent protein derived from Aequorea victoria. EGFP has F64L and S65T mutations which result in brighter green fluorescence.
产品形式	Liquid
产品应用	ELISA, WB, IP, IFA, FACS
反应性	Recognizes EGFP
克隆性	Monoclonal, clone 42-2
免疫原	Recombinant EGFP protein
浓度	1 mg/mL
标记	Unconjugated
产品组成	Phosphate Buffered Saline containing 0.03% ProClin 300, pH 7.4.
亚型	IgG2a

产品使用

Recommended dilution factor:

ELISA: 1:5000-20000

WB: 1:1000-5000

IP: 1:500-2000

IFA: 1:500-1000

FACS: 1:500-1000

Note: Working dilution for specific application should be determined by the investigator to obtain the best conditions.

产品储存/运输

储存

This product is stable after storage at: 2-8°C for 2 weeks under sterile conditions from date of receipt. -20°C or -80°C for 12 months under sterile conditions from date of receipt. Avoid repeated freeze/thaw cycles. Suggestion: Divide antibody into several vials. Keep only vials for usage at 2-8°C.

运输

The product is shipped with polar packs. Upon receipt, store it immediately at -20°C or lower for long term storage.

For Research Use or Further Manufacturing Only

Chamot Biotechnology(Shanghai) Co., Ltd. www.chamot-bio.com

Tel: 021-51880030 Mail: info@chamot-bio.com QQ: 864920491