

CHAMOT

**Transwell 小室 (基质胶包被) / Transwell Insert,
24-well plates with 8 μ m pores, PET, Sterile,
Matrix Gel-Coated**

CM003-12AR



CHAMOT

乔默® 生物

Specialize In Cytokines

CONTENT

1 产品简介

2 产品储存/运输

Transwell 小室（基质胶包被） / Transwell Insert, 24-well plates with 8 μ m pores, PET, Sterile, Matrix Gel-Coated

产品编号	CM003-12AR
规格	12 组/盒

产品简介

产品储存/运输

For Research Use or Further Manufacturing Only

Chamot Biotechnology(Shanghai) Co., Ltd. www.chamot-bio.com

Tel: 021-51880030 Mail: info@chamot-bio.com QQ: 864920491