CHAMOT

anti-Mouse CD3 Purified, 17A2 Clone, Functional assay

CM001-500MFA





トスツと

1 产品简介

2 产品储存/运输

CHAMOT 乔默[®]生物 Specialize In Cytokines

Product Datasheet

anti-Mouse CD3 Purified, 17A2 Clone, Functional assay

产品编号	CM001-500MFA
规格	500 μg

	And		
ᅳ	K	14	
10 10	1115	ו שו	

背景描述

Cluster of differentiation 3 (CD3) is one component of T cell receptor (TCR) complex that plays an important role in T cell development. CD3 cantains five different

that plays an important role in T cell development. CD3 cantains five different polypeptide chains: CD3-gamma, CD3-delta, CD3-epsilon, and zeta chain. The anti-CD3 antibody is specifically against the CD3-epsilon chain which is associated with

the antigen recognition and signal transduction.

别称 T-cell surface glycoprotein CD3

内毒素检测 <0.01 EU per 1 µg of the protein by the LAL method.

产品形式 Liquid

FC, FA, IHC, IP

产品应用 Note: It is recommended that the reagent be titrated for optimal performance for

each application.

反应性 Mouse

克隆性 Monoclonal, clone 17A2

浓度 2 mg/mL

纯化方式 Affinity Chromatography

标记 Unconjugated

产品组成 Phosphate-buffered aqueous solution, pH7.2, sodium azide-free

亚型 Rat IgG2b, kappa



The product should be stored undiluted at 2-8°C and protected from prolonged exposure to light. Do not freeze.



Blue Ice

For Research Use or Further Manufacturing Only

Chamot Biotechnology(Shanghai) Co., Ltd. www.chamot-bio.com

Tel: 021-51880030 Mail: info@chamot-bio.com QQ: 864920491