

# CHAMOT

## Donkey anti-Goat IgG (H+L)-HRP

CM005-A2HS

CM005-0.1A2H

CM005-0.5A2H

CM005-1A2H

CM005-5A2H

CM005-10A2H



CHAMOT

乔默®生物

Specialize In Cytokines



# CONTENT

1 产品简介

2 产品使用

3 产品储存/运输

## Donkey anti-Goat IgG (H+L)-HRP

产品编号	CM005-A2HS	CM005-0.1A2H	CM005-0.5A2H	CM005-1A2H	CM005-5A2H	CM005-10A2H
规格	10 μL	100 μL	500 μL	1 mL	5 mL	10 mL

### 产品简介

<b>产品描述</b>	Donkey anti-Goat IgG (H+L) Antibody (HRP) is conjugated with high-activity horseradish peroxidase (HRP).
<b>产品形式</b>	Liquid
<b>产品应用</b>	WB,ELISA,IHC
<b>反应性</b>	Goat
<b>克隆性</b>	Polyclonal
<b>免疫原</b>	Purified goat IgG
<b>浓度</b>	Lot-dependent, please refer to product label.
<b>标记</b>	HRP
<b>产品组成</b>	Phosphate Buffered Saline containing 1% BSA, 50% glycerol and 0.01%thimerosal, pH 7.4.
<b>亚型</b>	IgG

### 产品使用

Recommended dilution factor:

WB/ELISA: 1:5000-25000

IHC: 1:500-2500

Note: Working dilution for specific application should be determined by the investigator to obtain the best conditions.

## 产品储存/运输

This product is stable after storage at:

### 储存

2-8°C for 2 weeks under sterile conditions from date of receipt.  
-20°C or -80°C for 12 months under sterile conditions from date of receipt.

Avoid repeated freeze/thaw cycles and protect from light.

**Suggestion:** Divide antibody into several vials. Keep only vials for usage at 2-8°C.

### 运输

The product is shipped with polar packs. Upon receipt, store it immediately at -20°C or lower for long term storage.

For Research Use or Further Manufacturing Only

Chamot Biotechnology(Shanghai) Co., Ltd. [www.chamot-bio.com](http://www.chamot-bio.com)

Tel: 021-51880030 Mail: [info@chamot-bio.com](mailto:info@chamot-bio.com) QQ: 864920491