

# CHAMOT

## Donkey anti-Goat IgG (H+L)-HRP

CM005-A2HS

CM005-0.1A2H

CM005-0.5A2H

CM005-1A2H

CM005-5A2H

CM005-10A2H



CHAMOT

乔默®生物

Specialize In Cytokines



# CONTENT

1 产品简介

2 产品使用

3 产品储存/运输

## Donkey anti-Goat IgG (H+L)-HRP

|      |            |              |              |            |            |             |
|------|------------|--------------|--------------|------------|------------|-------------|
| 产品编号 | CM005-A2HS | CM005-0.1A2H | CM005-0.5A2H | CM005-1A2H | CM005-5A2H | CM005-10A2H |
| 规格   | 10 µL      | 100 µL       | 500 µL       | 1 mL       | 5 mL       | 10 mL       |

### 产品简介

#### 产品描述

Donkey anti-Goat IgG (H+L) Antibody (HRP) is conjugated with high-activity horseradish peroxidase (HRP).

#### 产品形式

Liquid

#### 产品应用

WB,ELISA,IHC

#### 反应性

Goat

#### 克隆性

Polyclonal

#### 免疫原

Purified goat IgG

#### 浓度

Lot-dependent, please refer to product label.

#### 标记

HRP

#### 产品组成

Phosphate Buffered Saline containing 1% BSA, 50% glycerol and 0.01%thimerosal, pH 7.4.

#### 亚型

IgG

### 产品使用

Recommended dilution factor:

WB/ELISA: 1:5000-25000

IHC: 1:500-2500

Note: Working dilution for specific application should be determined by the investigator to obtain the best conditions.

## 产品储存/运输

|    |                                                                                                                                                                                                                                                                                                                                                                            |
|----|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| 储存 | <p>This product is stable after storage at:</p> <p>2-8°C for 2 weeks under sterile conditions from date of receipt.<br/>-20°C or -80°C for 12 months under sterile conditions from date of receipt.</p> <p>Avoid repeated freeze/thaw cycles and protect from light.</p> <p><b>Suggestion:</b> Divide antibody into several vials. Keep only vials for usage at 2-8°C.</p> |
| 运输 | <p>The product is shipped with polar packs. Upon receipt, store it immediately at -20°C or lower for long term storage.</p>                                                                                                                                                                                                                                                |

For Research Use or Further Manufacturing Only

Chamot Biotechnology(Shanghai) Co., Ltd. [www.chamot-bio.com](http://www.chamot-bio.com)

Tel: 021-51880030    Mail: [info@chamot-bio.com](mailto:info@chamot-bio.com)    QQ: 864920491