

CHAMOT

Donkey anti-Goat IgG (H+L)-HRP

CM005-A2HS

CM005-0.1A2H

CM005-0.5A2H

CM005-1A2H

CM005-5A2H

CM005-10A2H



CHAMOT

乔默®生物

Specialize In Cytokines



CONTENT

1 产品简介

2 产品使用

3 产品储存/运输

Donkey anti-Goat IgG (H+L)-HRP

产品编号	CM005-A2HS	CM005-0.1A2H	CM005-0.5A2H	CM005-1A2H	CM005-5A2H	CM005-10A2H
规格	10 μ L	100 μ L	500 μ L	1 mL	5 mL	10 mL

产品简介

产品描述	Donkey anti-Goat IgG (H+L) Antibody (HRP) is conjugated with high-activity horseradish peroxidase (HRP).
产品形式	Liquid
产品应用	WB,ELISA,IHC
反应性	Goat
克隆性	Polyclonal
免疫原	Purified goat IgG
浓度	Lot-dependent, please refer to product label.
标记	HRP
产品组成	Phosphate Buffered Saline containing 1% BSA, 50% glycerol and 0.01%thimerosal, pH 7.4.
亚型	IgG

产品使用

Recommended dilution factor:

WB/ELISA: 1:5000-25000

IHC: 1:500-2500

Note: Working dilution for specific application should be determined by the investigator to obtain the best conditions.

产品储存/运输

This product is stable after storage at:

储存

2-8°C for 2 weeks under sterile conditions from date of receipt.
-20°C or -80°C for 12 months under sterile conditions from date of receipt.

Avoid repeated freeze/thaw cycles and protect from light.

Suggestion: Divide antibody into several vials. Keep only vials for usage at 2-8°C.

运输

The product is shipped with polar packs. Upon receipt, store it immediately at -20°C or lower for long term storage.

For Research Use or Further Manufacturing Only

Chamot Biotechnology(Shanghai) Co., Ltd. www.chamot-bio.com

Tel: 021-51880030 Mail: info@chamot-bio.com QQ: 864920491