

CHAMOT

Mouse anti-HA tag mAb

CM001-A1S

CM001-0.1A1

CM001-0.5A1

CM001-1A1

CM001-5A1

CM001-10A1



CHAMOT

乔默®生物

Specialize In Cytokines



CONTENT

1 产品简介

2 产品使用

3 产品储存/运输

Mouse anti-HA tag mAb

产品编号	CM001-A1S	CM001-0.1A1	CM001-0.5A1	CM001-1A1	CM001-5A1	CM001-10A1
规格	10 µg	100 µg	500 µg	1 mg	5 mg	10 mg

产品简介

产品描述

Mouse anti-HA tag mAb is specific to YPYDVDPDYA (HA) tag which can be applied to several immunoassays including ELISA, WB, IP, IFA and FACS. This antibody has been identified in several sources (bacteria, yeast and mammalian cells) of HA tag fusion protein.

产品形式

Liquid

产品应用

ELISA, WB, IP, IFA, FACS

反应性

Recognizes HA tag on either the N- or C-terminus of recombinant fusion protein.

克隆性

Monoclonal, clone 98-106

免疫原

Peptide: YPYDVDPDYA (HA)

浓度

1 mg/mL

标记

Unconjugated

产品组成

Phosphate Buffered Saline containing 50% glycerol and 0.02% sodium azide.

产品使用

Recommended dilution factor:

ELISA: 1:5000-20000

WB: 1:1000-5000

IP: 1:500-2000

IFA: 1:500-1000

FACS: 1:500-1000

Note: Working dilution for specific application should be determined by the investigator to obtain the best conditions.

产品储存/运输

储存	This product is stable after storage at: 2-8°C for 2 weeks under sterile conditions from date of receipt. -20°C or -80°C for 12 months under sterile conditions from date of receipt. Avoid repeated freeze/thaw cycles. Suggestion: Divide antibody into several vials. Keep only vials for usage at 2-8°C.
运输	The product is shipped with polar packs. Upon receipt, store it immediately at -20°C or lower for long term storage.

For Research Use or Further Manufacturing Only

Chamot Biotechnology(Shanghai) Co., Ltd. www.chamot-bio.com

Tel: 021-51880030 Mail: info@chamot-bio.com QQ: 864920491