

CHAMOT

**流感病毒A抗体 / Anti-Influenza A virus NP
Antibody [clone IA02]**

CM005-100AI



CHAMOT

乔默®生物

Specialize In Cytokines



CONTENT

1 产品简介

2 产品使用

3 产品储存/运输

流感病毒A抗体 / Anti-Influenza A virus NP Antibody [clone IA02]

产品编号	CM005-100AI
规格	100 µg

产品简介

产品描述 Anti-Influenza A virus NP Antibody only recognizes nucleocapsid protein (NP) of influenza A viruses but not influenza B viruses. Influenza A viruses are RNA viruses, and their subtypes are labeled according to an H number (for the type of hemagglutinin) and an N number (for the type of neuraminidase).

产品形式 Liquid

产品应用 ELISA

反应性 Influenza A virus Nucleocapsid protein

克隆性 Monoclonal, clone IA02

免疫原 Recombinant Influenza A virus nucleocapsid protein

浓度 1 mg/mL

标记 Unconjugated

产品组成 Phosphate Buffered Saline containing 0.03% ProClin 300, pH 7.4.

亚型 IgG1, IgG2a

产品使用

Recognize Influenza A virus NP in Lateral Flow, when clone IA02 antibody was paired with Anti-Influenza A virus NP Antibody [clone IA01].

Recommended dilution factor:

ELISA: 1:5000-20000

产品储存/运输

储存	The product is shipped with polar packs. Upon receipt, store it immediately at -20°C or lower for long term storage.
运输	This product is stable after storage at: 2-8°C for 2 weeks under sterile conditions from date of receipt. -20°C or -80°C for 12 months under sterile conditions from date of receipt. Avoid repeated freeze/thaw cycles. Suggestion: Divide antibody into several vials. Keep only vials for usage at 2-8°C.

For Research Use or Further Manufacturing Only

Chamot Biotechnology(Shanghai) Co., Ltd. www.chamot-bio.com

Tel: 021-51880030 Mail: info@chamot-bio.com QQ: 864920491