

CHAMOT

Goat anti-Rabbit IgG-iFluor 594

CM006-A2FS

CM006-0.1A2F

CM006-0.5A2F

CM006-1A2F

CM006-5A2F

CM006-10A2F



CHAMOT

乔默®生物

Specialize In Cytokines



CONTENT

1 产品简介

2 产品使用

3 产品储存/运输

4 文献引用

Goat anti-Rabbit IgG-iFluor 594

产品编号	CM006-A2FS	CM006-0.1A2F	CM006-0.5A2F	CM006-1A2F	CM006-5A2F	CM006-10A2F
规格	10 µL	100 µL	500 µL	1 mL	5 mL	10 mL

产品简介

产品描述 iFluor 594 dye is a bright, red-fluorescent dye that can be excited by the 561 nm or 594 nm laser. Leadgene Goat Anti-Rabbit IgG (H+L) Antibody (iFluor 594) is conjugated with high-activity iFluor 594.

产品形式 Liquid

产品应用 IF,FACS

反应性 Rabbit (Minor cross-reactivity to human in higher concentration)

克隆性 Polyclonal

免疫原 Purified rabbit IgG

浓度 0.5 mg/mL

标记 iFluor 594

产品组成 Phosphate Buffered Saline containing 1% BSA, 0.03% ProClin 300, pH 7.4.

亚型 IgG

产品使用

Recommended dilution factor:

IF: 1:200-1000

FACS:1:200-1000

Note: Working dilution for specific application should be determined by the investigator to obtain the best conditions.

产品储存/运输

储存

This product is stable after storage at: 2-8°C for 6 months under sterile conditions from date of receipt. Do not freeze and protect from light.

运输

The product is shipped with polar packs. Upon receipt, store it immediately at 2-8°C for long term storage.

文献引用

- Antagonizing apolipoprotein J chaperone promotes proteasomal degradation of mTOR and relieves hepatic lipid deposition. Hepatology. 2023 Oct 1;78(4):1182-1199. (IF 12.9)

For Research Use or Further Manufacturing Only

Chamot Biotechnology(Shanghai) Co., Ltd. www.chamot-bio.com

Tel: 021-51880030 Mail: info@chamot-bio.com QQ: 864920491