

CHAMOT

Mouse anti-mCherry mAb

CM008-A1S

CM008-0.1A1

CM008-0.5A1

CM008-1A1

CM008-5A1

CM008-10A1



CHAMOT

乔默®生物

Specialize In Cytokines



CONTENT

1 产品简介

2 产品使用

3 产品储存/运输

Mouse anti-mCherry mAb

产品编号	CM008-A1S	CM008-0.1A1	CM008-0.5A1	CM008-1A1	CM008-5A1	CM008-10A1
规格	10 µg	100 µg	500 µg	1 mg	5 mg	10 mg

产品简介

产品描述 The mCherry protein is a basic red fluorescent protein derived from *Discosoma* sp. Leadgene mouse anti-mCherry mAb can recognize native and denature mCherry protein but not GFP.

产品形式 Liquid

产品应用 ELISA, WB, IP, IFA, FACS

反应性 Recognizes mCherry but not GFP.

克隆性 Monoclonal, clone 45-5

免疫原 Recombinant mCherry protein

浓度 1 mg/mL

标记 Unconjugated

产品组成 Phosphate Buffered Saline containing 0.03% ProClin 300, pH 7.4.

亚型 IgG2a

产品使用

Recommended dilution factor:

ELISA: 1:5000-20000

WB: 1:1000-5000

IP: 1:500-2000

IFA: 1:500-1000

FACS: 1:500-1000

Note: Working dilution for specific application should be determined by the investigator to obtain the best conditions.

产品储存/运输

储存	This product is stable after storage at: 2-8°C for 2 weeks under sterile conditions from date of receipt. -20°C or -80°C for 12 months under sterile conditions from date of receipt. Avoid repeated freeze/thaw cycles. Suggestion: Divide antibody into several vials. Keep only vials for usage at 2-8°C.
运输	The product is shipped with polar packs. Upon receipt, store it immediately at -20°C or lower for long term storage.

For Research Use or Further Manufacturing Only

Chamot Biotechnology(Shanghai) Co., Ltd. www.chamot-bio.com

Tel: 021-51880030 Mail: info@chamot-bio.com QQ: 864920491