

CHAMOT

新冠病毒抗体 / COVID-19 Human Control Antibody Set (NP)(IgG/IgA/IgM)

CM013-0.1A1CoV



CHAMOT

乔默®生物

Specialize In Cytokines



CONTENT

1 产品简介

2 产品使用

3 产品储存/运输

新冠病毒抗体 / COVID-19 Human Control Antibody Set (NP) (IgG/IgA/IgM)

产品编号	CM013-0.1A1CoV
规格	100 µg

产品简介

产品描述 Human anti-SARS-CoV & CoV-2 NP antibody recognize human SARS Coronavirus Nucleocapsid protein (NP). Coronavirus NP localize to the cytoplasm and the nucleolus in both Virus-like particle (VLP) infected primary cells and in cells transfected with NP plasmid. NP has abundant expression in coronavirus and is a highly immunogenic phosphoprotein. NP is also conserved in sequence. Due to these characteristic above, the NP is an ideal marker for diagnosis.

产品形式 Liquid

产品应用 ELISA, WB, IFA, FACS

反应性 SARS-CoV & CoV-2 Nucleocapsid protein

克隆性 Recombinant Monoclonal

浓度 1 mg/mL

标记 Unconjugated

产品组成 Phosphate Buffered Saline, pH 7.4.

亚型 IgG, IgA, IgM

产品使用

Recommended dilution factor:

ELISA: 1:5000-20000

WB: Assay dependent

IFA: Assay dependent

FACS: Assay dependent

Note: Working dilution for specific application should be determined by the investigator to obtain the best conditions.

产品储存/运输

储存

The product is shipped with polar packs. Upon receipt, store it immediately at -20°C or lower for long term storage.

运输

This product is stable after storage at: 2-8°C for 2 weeks under sterile conditions from date of receipt. -20°C or -80°C for 12 months under sterile conditions from date of receipt. Avoid repeated freeze/thaw cycles. Suggestion: Divide antibody into several vials. Keep only vials for usage at 2-8°C.

For Research Use or Further Manufacturing Only

Chamot Biotechnology(Shanghai) Co., Ltd. www.chamot-bio.com

Tel: 021-51880030 Mail: info@chamot-bio.com QQ: 864920491