

# CHAMOT

**新冠病毒蛋白 / SARS-CoV-2 Nucleocapsid Protein  
(IHU B.1.640.2 Variant), Tag Free, E. coli**

CM021-100PCoV



CHAMOT

乔默®生物

*Specialize In Cytokines*



# CONTENT

1 产品简介

2 产品使用

3 产品储存/运输

## 新冠病毒蛋白 / SARS-CoV-2 Nucleocapsid Protein (IHU B.1.640.2 Variant), Tag Free, E. coli

|      |               |
|------|---------------|
| 产品编号 | CM021-100PCoV |
| 规格   | 100 µg        |

### 产品简介

#### 背景描述

SARS-CoV-2 is a kind of coronavirus which full name is severe acute respiratory syndrome coronavirus 2. SARS-CoV-2 is contagious that causes the respiratory diseases and lung diseases which make difficulty breathing. SARS-CoV-2 do the spillover event in 2019 because it has genetic diversity. SARS-CoV-2 is composed by four subunits (spike, envelope, membrane and nucleocapsid proteins). Its RNA genome is encapsulated with nucleocapsid protein. The viral envelope is comprised of spike, envelope and membrane protein. SARS-CoV-2 has high affinity to ACE2, which is highly expression in intestines, kidney, and heart.

#### 别称

Nucleoprotein, Nucleocapsid protein, NC Protein N

#### 蛋白编码

QSW58852.1

#### 分子量

The protein has a calculated MW of 45.68 kDa.

#### 表达系统

Escherichia coli

#### 纯度

>90% as determined by SDS-PAGE analysis.

#### 蛋白序列

Met1-Ala419

#### 产品形式

Lyophilized from a 0.2 µm filtered solution of PBS, pH 7.4.

### 产品使用

Centrifuge at 3000 rpm for 5 mins before opening. It is recommended to reconstitute the lyophilized protein in sterile H<sub>2</sub>O to a concentration not less than 100 µg/mL and incubate the stock solution at RT for at least 20 min to ensure sufficient re-dissolved.

## 产品储存/运输

### 储存

Lyophilized protein should be stored at -20°C for 1 year. Upon reconstitution, store at 2°C to 8°C for up to 1 week. Further dilute in a buffer containing a carrier protein or stabilizer (e.g. 0.1% BSA, 10%FBS, 5%HSA or 5% trehalose solution), protein aliquots should be stored at -20°C or -80°C for 3-6 months. Avoid repeated freeze/thaw cycles.

### 运输

Blue Ice

For Research Use or Further Manufacturing Only

Chamot Biotechnology(Shanghai) Co., Ltd. [www.chamot-bio.com](http://www.chamot-bio.com)

Tel: 021-51880030    Mail: [info@chamot-bio.com](mailto:info@chamot-bio.com)    QQ: 864920491