

CHAMOT

**重组小鼠白细胞介素33 /Recombinant IL-33,
Mouse, Animal-Free**

CM030-5MP

CM030-20MP

CM030-100MP

CM030-500MP

CM030-1000MP



CHAMOT

乔默®生物

Specialize In Cytokines



CONTENT

1 产品简介

2 产品使用

3 产品储存/运输

4 实验数据展示

重组小鼠白细胞介素33 /Recombinant IL-33, Mouse, Animal-Free

| 产品编号 | CM030-5MP | CM030-20MP | CM030-100MP | CM030-500MP | CM030-1000MP |
|------|-----------|------------|-------------|-------------|--------------|
| 规格 | 5 µg | 20 µg | 100 µg | 500 µg | 1 mg |

产品简介

背景描述

Interleukin 33 (IL-33) is a 17.65 kDa cytokine with 159 amino acid residues. IL-33 is a member of the IL-1 family and secreted from various cell types, such as mast cells, macrophages, endothelial cells, and epithelial cells. IL-33 activates NF-κB and MAPK signaling pathways that stimulate the secretion of type 2 cytokines like IL-5 and IL-13 when IL-33 binds to the ST2 receptor primarily expressed in mast cells and Th2 cells. In addition to Th type 2 responses, IL-33 participates in maintaining barrier tissue defense.

别称

interleukin 33, IL-1 F11, NF-HEV, 9230117N10Rik, IL33

蛋白编码

Q8BVZ5

分子量

The protein has a calculated MW of 18.51 kDa. The protein migrates as 17-25 kDa under reducing condition (SDS-PAGE analysis).

表达系统

Escherichia coli

纯度

>98% as determined by SDS-PAGE.

生物活性

Measure by its ability to induce proliferation in D10.G4.1 cells. The ED₅₀ for this effect is <40 pg/mL. The specific activity of recombinant mouse IL-33 is > 2 x 10⁷ IU/mg.

内毒素检测

<0.1 EU per 1 µg of the protein by the LAL method.

蛋白序列

MSIQGTSLLTQSPASLSTYNDQSVSFVLENGCYVINVDDSGKDQEQDQVLLRYYESPCPASQS
GDGVDGKKLMVNMSPDKTDIWLHANDKDYSVELQRGDVSPPEQAFFVLHKKSSDFVSFEC
KNLPGTYIGVKDNQLALVEEKDESCNNIMFKLSKI with polyhistidine tag at the C-terminus.

产品形式

The protein was lyophilized from a 0.2 µm filtered solution containing 1X PBS, pH 7.4. If you have any concerns or special requirements, please confirm with us.

产品应用

Cell Culture

产品使用

1. Before opening, centrifuge at 3000 rpm for 5 mins.

2.Initial Reconstitution

- Reconstitute the Lyophilized Protein in sterile H₂O to a concentration of 100-200 µg/mL.
- Then, incubate it at room temperature for at least 20 mins to ensure sufficient dissolution.

Do Not Vortex! Vigorous shaking may impair the biological activity of the protein.

- Store at 2°C to 8°C for up to 1 week.

3.Extended Storage

- After the initial reconstitution, further dilute the reconstituted protein in a buffer containing a carrier protein or stabilizer (e.g., 0.1% BSA, 10% FBS, 5% HSA, or 5% trehalose solution). The final concentration is not less than 10 µg/ml.

- Prepare aliquots (≥20 µl).
- Store at -20°C or -80°C for 3 to 6 months.

Avoid repeated freeze-thaw cycles.

产品储存/运输

Lyophilized Protein: Store at -20°C for 1 year.

储存


Initial Reconstitution: Store at 2°C to 8°C for up to 1 week.

Extended Storage: Store at -20°C or -80°C for 3 to 6 months with a carrier protein or stabilizer.

运输

Blue Ice

实验数据展示



For Research Use or Further Manufacturing Only

Chamot Biotechnology(Shanghai) Co., Ltd. www.chamot-bio.com

Tel: 021-51880030 Mail: info@chamot-bio.com QQ: 864920491